

MPE ASR


Automatic Condensing Waste Discharge System

Versatile Condensation Receiving Module

- For mixed solvent with high and low boiling point in evaporation, automatic condensing waste discharge system can discharge low-boiling solvent first, to avoid interference with the subsequent samples.

Precise Pressure Control

- Available in various built-in pressure-decrease modes: manual, programmed, frequency-conversion.
- Flexible to different requirements: able to support for 20-step gradient control, with optional frequency-conversion pump, able to control gradient change of pressure smoothly.



RayKol Group

Add: 5-6F, No.176 Xinfeng Road, HuizhiZone, Torch High-tech Zone, Xiamen, China

Tel: +86 0592 5800190

Mail: info@raykolgroup.com

www.raykolgroup.com